

**FT458-FO BodyGuard Wall Padding**

1. **PADDING DESIGN** – Each wall pad section shall be constructed using foam with vinyl covering and velcro attachment. Foam shall be 2” thick wrapped in a heavy duty 14 oz. vinyl. A 2” wide velcro strip shall be sewn horizontally to the top and bottom of the pad. A mating piece of velcro shall be supplied with an adhesive back for mounting to wall surface. Additional hardware may be recommended for securing velcro to wall where adhesive back may be insufficient (wall fastening hardware not provided). Overall pad size is 24” wide by 96” tall. Available in 16 colors.